

DJL-8/12/16/24/32 Electronic Tablets Counting Machine



Executive Standard:

- (GB/T 29020-2012)) (Our company is the main unit)
- (JB/T20019-2014)) (Our company is the revision unit)

High adaptability:

- A pill with a minimum diameter 2.5mm.
- Apply various type tablets, including special-shaped tablets.
- Apply soft&hard capsules, transparent and translucent capsules.
- Vibration blanking: 8 channels vibration uniform blanking, unique patent plate loading mechanism blanking, uniform and stable discharging, no damage to drugs.
- Anti-high dust: the company's original anti-high dust photoelectric sensing technology can work stably under the condition of high dust.
- Accurate counting: photoelectric sensor automatic counting, bottling error less than the national standard.
- Leakage nozzle with special structure can prevent clogging, and the bottle with small mouth can be filled quickly (optional).
- High degree of intelligence, with no bottles, no number, abnormal bottles or a certain size of the broken grain automatic elimination, fault self-check and other detection and alarm control functions.
- Supporting functions: the lack of bottle photoelectric sensor automatic detection, automatic shutdown.

- Simple operation: intelligent design, Chinese or English menu, various operating parameters according to the need to set, can store 10 groups of setting parameters, easy to call.
- Convenient maintenance: operators only need simple training to operate, disassembly, cleaning, replacing parts is simple and convenient, without tools.
- Sealed and dustproof: Equipped with dust box to reduce dust pollution (optional).

Parameter Information

Model	DJL-8	DJL-12	DJL-16
Capacity	10-45 Bottles per minute	10-45 Bottles per minute	20-80 Bottles per minute
Bottle quantity setting range	2-9999		
Voltage and power	220V/50HZ 0.7kw	220V/50HZ 0.8kw	220V/50HZ 1.2kw
Size	1200×1280×1560mm	1400×1400×1560mm	1600×1400×1560mm
Filling error	≤ 2‰	≤ 2‰	≤ 2‰
Main material of equipment	SS316L(The part in contact with the drug) +SS304	SS316L(The part in contact with the drug) +SS304	SS316L(The part in contact with the drug) +SS304
Control system	MCU	MCU	MCU
The compressed air	0.6 Mpa	0.6 Mpa	0.6 Mpa
Weight	180kg	270kg	410kg

Model	DJL-24	DJL-32	DJL-48
Capacity	20-90 Bottles per minute	40-120 Bottles per minute	40-150 Bottles per minute
Bottle quantity setting range	2-9999		
Voltage and power	220V/50HZ 1.5kw	220V/50HZ 2.2kw	220V/50HZ 2.5kw
Size	2200×1400×1680mm	2200×1400×1680mm	2200×1400×1680mm
Filling error	≤ 2‰	≤ 2‰	≤ 2‰
Main material of equipment	SS316L(The part in contact with the drug) +SS304	SS316L(The part in contact with the drug) +SS304	SS316L(The part in contact with the drug) +SS304
Control system	MCU	MCU	MCU
The compressed air	0.6 Mpa	0.6 Mpa	0.6 Mpa
Weight	460kg	650kg	750kg