

MP Effervescent Tablets Counting & Filling Machine



Specially designed for the packaging of straight tube for effervescent tablets, the tablets can be packed in a straight tube with a set amount of the numbers.

- Vibration blanking: multi-stage vibration uniform blanking, special mechanism blanking, uniform and stable discharging, no damage to drugs
- Main function: put the tablets into the straight bottle according to the set quantity.
- Detection function: with the missing bottle alarm does not fill, broken particles automatically remove; A number of detection and alarm control functions, such as fixed-point elimination and fault self-check, when the loading quantity does not meet the requirements.
- Anti-high dust: adopt unique anti-high dust photoelectric sensor.
- Special structure of the bottle system, mold bottle.
- Supporting functions: the lack of bottle photoelectric sensor automatic detection, automatic shutdown.
- Simple operation: intelligent design, man-machine dialogue, Chinese menu, various operating parameters according to the need to set, can store 10 groups of setting parameters, in order to call.
- Convenient maintenance: operators only need simple training to operate, disassembly, cleaning, replacing parts is simple and convenient, without tools.
- Sealed and dust-proof: Equipped with dust box to reduce dust pollution (optional).

Technical Parameters

Model	MP
Tablets specifications	Each group is 10-16 in diameter and 18-25 in diameter
Capacity	40-80 BPM
Voltage and power	220V/50HZ 2.0 kw
The compressed air	0.6 Mpa
Size	2620×2180×1950mm
Weight	600kg