

## WJ-1 Dropping Pills Counting & Filling Machine



- It is mainly used to pack the 2.5mm-5mm dropping pills, especially for special shaped bottles. Bottle unscrambling unit,
- counting and filling unit, capping unit processes can be done in one machine. Solve the problem of small workshop area, save labor, only need 1 person operation.

### Technical Parameters

Model	WJ-1
Capacity	Take 30 capsules/bottle as an example, 30-50 bottles/minute, 13 500-22 500 bottles/shift
Voltage and power	220V/50HZ 0.8-1.2kw
Size	1500×1500×1500mm
Main material of equipment	SS316+SS304
The compressed air	0.6 Mpa
Weight	400-500kg