



PY-85 Drug Deblister Machine



PY-85 Drug Deblister Machine is applied to eliminate and recycle capsule, tablet, sugar-coated tablet and so on. To solve the plate finished product capsule, tablet, sugar pill, and the separation of aluminum-plastic plate, to reduce the labor strength and improve the working efficiency. Make the whole recycling process more standard, simplified, greatly improve production efficiency.

Technical features

- a) Simple structure, easy operation, does not need professional operation
- b) High efficiency and reliable operation
- c) Made by stainless steel is completely sealed
- d) Speed is adjustable, the infinite speed

Technical Parameter

| Model | PY-85 |
|---------------------------|----------------|
| Production Capacity | 65 pieces/min |
| Maximum adjustment size | 150*120mm |
| Maximum number of columns | 7 columns |
| Total Power | 90W |
| Power | 220V, 50HZ, 1A |
| Weigh | 78kg |
| Dimension | 520×420×1300mm |